


- 
- Flexible Hydrocolloid adhesive
 - Soft feeling noiseless pouch film
 - Extra-thin and dry touch nonwoven liner
 - Secure coupling system for two-piece


WE DO OSTOMY CARE



BAO-HEALTH MEDICAL INSTRUMENT CO., LTD.

An'Guard/ One Piece Pouch

For colostomy/Ileostomy stoma



Code No. : 1023TNF

One Piece Pouch (with clamp closure)

Another choice of daily care for col/Ile stoma.

- Skin friendly
- Good adhesive ability
- Good anti-leakage function
- Odor free and noiseless

Code No.: 1023TNF

Pouch size: 15X28cm Cutting range: 15-60mm

Pouch volume: 650ml Clamp closure

CureGuard/ One Piece Pouch

Classical pouch

Excellent feeling from Cureguard adhesive

The Cureguard adhesive is thin and elastic in nature which enables it to follow the contours and movements of the skin.

Meanwhile it acts as a skin barrier, prevent the output to get in contact with the skin and reducing skin irritation.

With a gentle tack and high absorb moisture rate, Cureguard assures long wear life. High security against leakage and complete skin protection.




Code No.: 1145001F

Pouch size: 15X28cm

Pouch volume: 650ml2.

Cutting range: 15-60mm

Magic closure



Code No.: 1117001F

Pouch size: 15X28cm

Pouch volume: 650ml2.

Cutting range: 15-70mm

Integral flex closure

XL/XXL pouch



For post-operation use when stoma swell,
or for loop-stoma daily care

Code No.: 10875001F

Pouch size: 16X30cm

Pouch volume: 750ml

Cutting range: 15-75mm

Code No.: 10910201F


Pouch size: 16X30cm

Pouch volume: 750ml

Cutting range: 15-102mm


CureGuard/ Two Piece System

Classical pouch



- Cureguard adhesive skin barrier.
- Flexible PE coupling system, anti-leakage design.
- Thicker Cureguard, prolong using time.
- Perfect anti-leakage coupling system, no leakage will occur even under high pressure(20kgs), and bending.

Code No.:	3015001	3016001	3017001
Cutting Range:	15-50mm	15-60mm	15-70mm
Plate size:	10cm x 10cm	10cm x 10cm	11.5cm x 11.5cm



Code No.: 2025008F

Pouch size: 15X28cm

Pouch volume: 650ml

Match to:3015001

Clamp closure


Code No.: 2027008F

Pouch size: 15X29cm

Pouch volume: 700ml

Match to:3017001

Clamp closure



Code No.: 2146008F

Pouch size: 15X28cm

Pouch volume: 650ml

Match to:3016001


Magic closure

NOTE: SKIN BARRIER AND POUCHES COULD BE SIZE MATCHED BY SAME COLOR



Code No.: 3064501

45mm



"ZERO" Pressure Skin barrier NEW FOR 2017

- Cureguard adhesive skin barrier.
- Flexible PE coupling system, **pressure-free**.
- Thicker Cureguard, prolong using time.
- Perfect anti-leakage coupling system, no leakage will occur even under high pressure(20kgs), and bending.


1.No pressure on belly when install pouch.

2.Rotatable design for coupling system

3.Softor feeling, more comfortable.

CureGuard/Urostomy bag

One Piece Pouch for Urostomy



Products Features:


- Cureguard hydrocolloid has very good anti-leakage function.
- Odor free and noiseless.
- Soft and dry touch of nonwoven liner.
- Elastic drainage easy to control to drain the urine, prevent dirty hands.
- Anti-flow back design.

Code No.: 4016001

Pouch size: 15X23cm Cutting range: 15-60mm

Pouch volume: 550ml

Two Piece Pouch for Urostomy



Code No.: 4035001

Pouch size: 14X23cm

Pouch volume: 550ml

Match to: 3015001

Code No.: 4036001


Pouch size: 14X23cm

Pouch volume: 550ml

Match to: 3016001

NOTE: SKIN BARRIER AND POUCHES COULD
BE SIZE MATCHED BY SAME COLOR

**PEOPLE IN 14 COUNTRIES
NOW USE BAO-HEALTH**



CE FDA



Certificates

Exclusive Distributor in India: