



Grade	Compression Intensity	Pressure in mmHg	Indication
I	Low	18-21	Feeling of tired & heavy leg, early venous insufficiency (weakness), mild varicose veins without pronounced edema, varicose veins in pregnancy, spider veins
II	Moderate	23-32	More serious symptoms, more marked varicose veins or moderate edema, after healing of minor ulcers, after scleriosis and surgery for varicose veins, DVT, in sufficiency of the valves in the deep veins
III	High	34-46	Chronic Venous Insufficiency, pronounced edema, atrophie blanche, dermatosclerosis, after resolution of severe ulcers
IV	Very High	49 & Higher	Lymphodema, disfiguring swelling






Measurement Chart for Thigh High Stockings							
	XS	S	M	L	XL	XXL	XXXL
A / Minimum Ankle, cm	15 - 17	18 - 20	21 - 23	24 - 26	27 - 29	30 - 32	33 - 36
B / Maximum Calf, cm	24 - 30	28 - 34	32 - 38	36 - 42	40 - 46	44 - 50	48 - 54
C / Mid Thigh, cm	37 - 50	42 - 57	47 - 64	54 - 71	61 - 78	68 - 85	75 - 92
T / Thigh Length cm	72 - 83	72 - 83	72 - 83	72 - 83	72 - 83	72 - 83	76 - 85



Measurement Chart for Knee High Stockings & Wound Care Socks							
	XS	S	M	L	XL	XXL	XXXL
A / Minimum Ankle, cm	15 - 17	18 - 20	21 - 23	24 - 26	27 - 29	30 - 32	33 - 36
B / Maximum Calf, cm	24 - 30	28 - 34	32 - 38	36 - 42	40 - 46	44 - 50	48 - 54
T / Knee Length, cm	39 - 44	39 - 44	39 - 44	42 - 48	42 - 48	42 - 48	44 - 52



Measurement Chart for Arm Sleeve						
	S	M	L	XL	XXL	
A / Wrist, cm	15 - 16	16 - 19	18 - 21	20 - 23	22 - 25	
B / Bicep, cm	19 - 25	22 - 29	25 - 32	29 - 36	33 - 40	
T / Arm Length	39 - 44	39 - 44	39 - 44	44 - 49	44 - 49	



## Your Local Distributor



Imported & Marketed in India by :  
**KANAV HEALTHCARE PVT. LTD.**  
 Sushant Lok-1, Gurgaon 122009 Haryana  
 +91 124 4273077  
 kanavhealthcare@gmail.com  
 matrisurg.com



AWARDED THE BEST BRAND OF COMPRESSION STOCKINGS

